


RASAL1 Antibody (R35254)

Catalog No.	Formulation	Size
R35254-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide $$	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	8437
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This RASAL1 antibody is available for research use only.



Description

Additional name(s) for this target protein: RAS protein activator like 1, RASAL

Application Notes

Optimal dilution of the RASAL1 antibody should be determined by the researcher.

Immunogen

Amino acids QERGKAALGPLGP were used as the immunogen for this RASAL1 antibody.

Storage Aliquot and store the RASAL1 antibody at -20oC.