


Ribosomal Protein L19 Antibody (R34738)

Catalog No.	Formulation	Size
R34738-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat, Dog
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	6143
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This Ribosomal Protein L19 antibody is available for research use only.



Description

Additional name(s) for this target protein: RPL19; L19

Application Notes

Optimal dilution of the Ribosomal Protein L19 antibody should be determined by the researcher.

Immunogen

Amino acids LRRYRESKKIDRH were used as the immunogen for this Ribosomal Protein L19 antibody. **Storage** Aliquot and store the Ribosomal Protein L19 antibody at -20oC. Ordering: Phone: 858.663.9055 | Fax: 1.267.821.0800 | Email: info@nsjbio.com $\label{lem:copyright @NSJ Bioreagents.} \ \ \text{All rights reserved}.$