


## RICTOR Antibody (F50244)

Catalog No.	Formulation	Size
F50244-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F50244-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

<b>Availability</b>	1-3 business days
<b>Species Reactivity</b>	Human, Mouse
<b>Format</b>	Purified
<b>Clonality</b>	Polyclonal (rabbit origin)
<b>Isotype</b>	Rabbit Ig
<b>Purity</b>	Purified
<b>UniProt</b>	Q6R327
<b>Applications</b>	Western Blot : 1:1000
<b>Limitations</b>	This RICTOR antibody is available for research use only.



### Description

RICTOR and MTOR (FRAP1) are components of a protein complex that integrates nutrient- and growth factor-derived signals to regulate cell growth.

### Application Notes

Titration of the RICTOR antibody may be required due to differences in protocols and secondary/substrate sensitivity.

### Immunogen

A portion of amino acids 805-835 from the human protein was used as the immunogen for this RICTOR antibody.

## Storage

Aliquot the RICTOR antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.