


Sox17 Antibody (F53214)

Catalog No.	Formulation	Size
F53214-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F53214-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Format	Antigen affinity purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	Q61473
Applications	Western Blot : 1:1000-2000
Limitations	This Sox17 antibody is available for research use only.



Western blot testing of Sox17 antibody at 1:2000 dilution. Lane 1: mouse F9 lysate; 2: human LNCaP lysate; 3: mouse testis lysate; Predicted molecular weight ~45 kDa.

Description

Sox17 acts as transcription regulator that binds target promoter DNA and bends the DNA.

Application Notes

Titration of the Sox17 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 350-400 (mouse) was used as immunogen for this Sox17 antibody.

Storage

Aliquot the Sox17 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.