


SRD5A1 Antibody (R33183)

| Catalog No. | Formulation | Size |
|--------------|---|--------|
| R33183-100UG | 0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide | 100 ug |

[Bulk quote request](#)

| | |
|--------------------|---|
| Availability | 1-3 business days |
| Species Reactivity | Human |
| Format | Antigen affinity purified |
| Clonality | Polyclonal (goat origin) |
| Isotype | Goat Ig |
| Purity | Antigen affinity |
| Gene ID | 6715 |
| Applications | Immunohistochemistry (FFPE) : 4-6ug/ml ELISA (peptide) LOD : 1:64000 |
| Limitations | This SRD5A1 antibody is available for research use only. |



IHC testing of FFPE human prostate with SRD5A1 antibody at 4ug/ml. HIER: steamed with pH9 Tris/EDTA buffer, HRP-staining.

Description

Additional name(s) for this target protein: Steroid-5-alpha-reductase, alpha polypeptide 1

Application Notes

Optimal dilution of the SRD5A1 antibody should be determined by the researcher.

Immunogen

Amino acids ATATGVAEERLLC were used as the immunogen for this SRD5A1 antibody.

Storage

Aliquot and store the SRD5A1 antibody at -20oC.