


Stap2 Antibody (R33749)

Catalog No.	Formulation	Size
R33749-100UG	$0.5\ \text{mg/ml}$ in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Mouse
Predicted Reactivity	Cow, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	106766
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This Stap2 antibody is available for research use only.



Description

Additional name(s) for this target protein: Signal transducing adaptor family member 2

Application Notes

Optimal dilution of the Stap2 antibody should be determined by the researcher.

Immunogen

Amino acids EEKLQKRRALEH were used as the immunogen for this Stap2 antibody. **Storage** Aliquot and store the Stap2 antibody at -20oC.