


## STAT5A Antibody (R33643)

Catalog No.	Formulation	Size
R33643-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	6776
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This STAT5A antibody is available for research use only.



Western blot testing of human A) K562, B) MOLT4 and C) mouse NIH3T3 lysate with STAT5A antibody at 0.3ug/ml. Predicted molecular weight: ~91kDa.

## Description

Additional name(s) for this target protein: Signal transducer and activator of transcription 5A

## Application Notes

Optimal dilution of the STAT5A antibody should be determined by the researcher.

1. This STAT5A antibody will not cross react with STAT5B.

## Immunogen

Amino acids DSRLSPPAGLFTSAR were used as the immunogen for this STAT5A antibody.

## Storage

Aliquot and store the STAT5A antibody at -20oC.