


STEAP4 Antibody (R33360)


| Catalog No. | Formulation | Size |
|--------------|---|--------|
| R33360-100UG | 0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide | 100 ug |

[Bulk quote request](#)

| | |
|----------------------|---|
| Availability | 1-3 business days |
| Species Reactivity | Human |
| Predicted Reactivity | Dog, Mouse, Rat |
| Format | Antigen affinity purified |
| Clonality | Polyclonal (goat origin) |
| Isotype | Goat Ig |
| Purity | Antigen affinity |
| Gene ID | 79689 |
| Applications | Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:4000 |
| Limitations | This STEAP4 antibody is available for research use only. |



Western blot testing of human adipose lysate with STEAP4 antibody at 0.3ug/ml.
Predicted molecular weight: ~52kDa.



Western blot of HEK293 lysate overexpressing STEAP4 probed with STEAP4 antibody (mock transfection in lane 1). Predicted molecular weight: ~52kDa.

Description

Additional name(s) for this target protein: STEAP family member 4; tumor necrosis factor, alpha-induced protein 9; TNFAIP9

Application Notes

Optimal dilution of the STEAP4 antibody should be determined by the researcher.

Immunogen

Amino acids CVDNTLTRIRQGWERN were used as the immunogen for this STEAP4 antibody.

Storage

Aliquot and store the STEAP4 antibody at -20oC.