

STK19 Antibody (F40147)

Catalog No.	Formulation	Size
F40147-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F40147-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Purified
UniProt	P49842
Applications	Western Blot : 1:1000
Limitations	This STK19 antibody is available for research use only.

130
72
56
36
28
17
11

Western blot analysis of STK19 antibody and HL-60 lysate.

Description

SSTK19/STK19 is a serine/threonine kinase which localizes predominantly to the nucleus. Its specific function is unknown; it is possible that phosphorylation of this protein is involved in transcriptional regulation.

Application Notes

Titration of the STK19 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 278-308 from the human protein was used as the immunogen for this STK19 antibody.

Storage

Aliquot the STK19 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.