


STK25 Antibody (F54066)

Catalog No.	Formulation	Size
F54066-0.2ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.2 ml
F54066-0.05ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.05 ml

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	O00506
Applications	Western Blot : 1:2000
Limitations	This STK25 antibody is available for research use only.



Western blot testing of STK25 antibody at 1:2000: Lane 1) human A431 cell lysate, 2) human HepG2 cell lysate and 3) human testis lysate. Predicted molecular weight ~48 kDa.

Description

Oxidant stress-activated serine/threonine kinase that may play a role in the response to environmental stress. Targets to the Golgi apparatus where it appears to regulate protein transport events, cell adhesion, and polarity complexes important for cell migration. [UniProt]

Application Notes

The stated application concentrations are suggested starting points. Titration of the STK25 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

Recombinant human STK25 was used as the immunogen for the STK25 antibody.

Storage

Aliquot the STK25 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.