

SUV39H2 Antibody (R33257)

Catalog No.	Formulation	Size
R33257-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Predicted Reactivity	Dog, Pig, Cow
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	79723
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This SUV39H2 antibody is available for research use only.



Western blot testing of rat kidney lysate with SUV39H2 antibody at 1ug/ml. Expected molecular weight: 40-53 kDa.

Description

Additional name(s) for this target protein: Suppressor of variegation 3-9 homolog 2

Application Notes

Optimal dilution of the SUV39H2 antibody should be determined by the researcher.

Immunogen

Amino acids CKCGAVTCRGYLN were used as the immunogen for this SUV39H2 antibody.

Storage

Aliquot and store the SUV39H2 antibody at -20oC.