


# Synaptotagmin-9 Antibody (R33794)

Catalog No.	Formulation	Size
R33794-100UG	0.5~mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

#### **Bulk quote request**

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Predicted Reactivity	Cow, Pig
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	143425
Applications	Western Blot : 0.5-2ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This Synaptotagmin-9 antibody is available for research use only.



## **Description**

Additional name(s) for this target protein: SYT9; SYTIX

### **Application Notes**

Optimal dilution of the Synaptotagmin-9 antibody should be determined by the researcher.

### **Immunogen**

Amino acids EHDSCQDFIYHLRD were used as the immunogen for this Synaptotagmin-9 antibody. **Storage** Aliquot and store the Synaptotagmin-9 antibody at -20oC. Ordering: Phone: 858.663.9055 | Fax: 1.267.821.0800 | Email: info@nsjbio.com  $\label{lem:copyright @NSJ Bioreagents.} \ \ \text{All rights reserved}.$