


SYNJ2BP Antibody (F53989)

Catalog No.	Formulation	Size
F53989-0.2ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.2 ml
F53989-0.05ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.05 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat
Format	Antigen affinity purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	P57105
Applications	Western Blot : 1:1000-2000
Limitations	This SYNJ2BP antibody is available for research use only.



Western blot testing of human PC-3 cell lysate with SYNJBP antibody at 1:2000.
Predicted molecular weight: 16 kDa.

Description

Additional name(s): Synaptojanin 2 binding protein

Application Notes

Titration of the SYNJ2BP antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 22-66 from the human protein was used as the immunogen for the SYNJ2BP antibody.

Storage

Aliquot the SYNJ2BP antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.