


# SYNJ2BP Antibody (F53989)

Catalog No.	Formulation	Size
F53989-0.2ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.2 ml
F53989-0.05ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.05 ml

## **Bulk quote request**

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	P57105
Applications	Western Blot : 1:1000-2000
Limitations	This SYNJ2BP antibody is available for research use only.



#### **Description**

Additional name(s): Synaptojanin 2 binding protein

# **Application Notes**

Titration of the SYNJ2BP antibody may be required due to differences in protocols and secondary/substrate sensitivity.

### **Immunogen**

A portion of amino acids 22-66 from the human protein was used as the immunogen for the SYNJ2BP antibody. **Storage** Aliquot the SYNJ2BP antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles. Ordering: Phone: 858.663.9055 | Fax: 1.267.821.0800 | Email: info@nsjbio.com Copyright  $\ensuremath{\texttt{@}}$  NSJ Bioreagents. All rights reserved.