

TCF4 Antibody (R34798)

| Catalog No. | Formulation | Size |
|--------------|---|--------|
| R34798-100UG | 0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide | 100 ug |

[Bulk quote request](#)

| | |
|----------------------|--|
| Availability | 1-3 business days |
| Species Reactivity | Human |
| Predicted Reactivity | Dog, Mouse |
| Format | Antigen affinity purified |
| Clonality | Polyclonal (goat origin) |
| Isotype | Goat Ig |
| Purity | Antigen affinity |
| Gene ID | 6925 |
| Applications | Western Blot : 0.5-1.5ug/ml ELISA (peptide) LOD : 1:16000 |
| Limitations | This TCF4 antibody is available for research use only. |



Western blot testing of human heart lysate with TCF4 antibody at 0.5ug/ml. Predicted molecular weight ~71kDa.

Description

Additional name(s) for this target protein: Transcription factor 4, E2-2, ITF2

Application Notes

Optimal dilution of the TCF4 antibody should be determined by the researcher.

1. This TCF4 antibody will detect isoforms a and b.

Immunogen

Amino acids NGGHPSPSRNYGDGT were used as the immunogen for this TCF4 antibody.

Storage

Aliquot and store the TCF4 antibody at -20oC.