


TDP1 Antibody (R35075)

Catalog No.	Formulation	Size
R35075-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	55775
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This TDP1 antibody is available for research use only.



Western blot testing of Jurkat nuclear lysate with TDP1 antibody at 2ug/ml. Predicted molecular weight: ~68kDa.

Description

Additional name(s) for this target protein: Tyrosyl-DNA phosphodiesterase 1

Application Notes

Optimal dilution of the TDP1 antibody should be determined by the researcher.

Immunogen

Amino acids PNDGTAQRTEHNG were used as the immunogen for this TDP1 antibody.

Storage

Aliquot and store the TDP1 antibody at -20oC.