


THEM2 Antibody [clone 169CT7.1.3] (F40235)

Catalog No.	Formulation	Size
F40235-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F40235-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Purified
Host	Mouse
Clonality	Monoclonal (mouse origin)
Isotype	Mouse IgG1, k
Clone Name	169CT7.1.3
Purity	Purified
UniProt	Q9NPJ3
Applications	Western Blot : 1:100-1:1000
Limitations	This THEM2 antibody is available for research use only.



Description

THEM2 belongs to the UPF0152 (paal) family. Its exact function is still undetermined.

Application Notes

Titration of the THEM2 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

This THEM2 antibody was produced from a mouse immunized with THEM2 recombinant protein.

Storage

Aliquot the THEM2 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.