


## Themis Antibody (R33761)

Catalog No.	Formulation	Size
R33761-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Mouse
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	210757
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This Themis antibody is available for research use only.



Western blot testing of mouse spleen lysate with Themis antibody at 0.3ug/ml. Predicted molecular weight: ~73kDa.

## Description

Additional name(s) for this target protein: Thymocyte-expressed molecule involved in selection

## Application Notes

Optimal dilution of the Themis antibody should be determined by the researcher.

## Immunogen

Amino acids EERTINLPKSLKSH were used as the immunogen for this Themis antibody.

## Storage

Aliquot and store the Themis antibody at -20oC.