


## TIP30 Antibody (R33127)

Catalog No.	Formulation	Size
R33127-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
UniProt	Q9BUP3
Gene ID	10553
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This TIP30 antibody is available for research use only.



Western blot testing of HeLa lysate with TIP30 antibody at 0.6ug/ml. Predicted molecular weight: ~27/15/30kDa (isoforms 1/2/3).

## Description

Additional name(s) for this target protein: HIV-1 Tat interactive protein 2, 30kDa; HTATIP2

## Application Notes

Optimal dilution of the TIP30 antibody should be determined by the researcher.

## Immunogen

Amino acids AETEALSKLREDFR were used as the immunogen for this TIP30 antibody. This sequence is common to isoforms 1-3.

## Storage

Aliquot and store the TIP30 antibody at -20oC.