


TMEM205 Antibody (R34979)

Catalog No.	Formulation	Size
R34979-100UG	$0.5\ \text{mg/ml}$ in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	374882
Applications	Western Blot: 0.03-0.1ug/ml ELISA (peptide) LOD: 1:16000
Limitations	This TMEM205 antibody is available for research use only.



Description

Additional name(s) for this target protein: Transmembrane protein 205

Application Notes

Optimal dilution of the TMEM205 antibody should be determined by the researcher.

Immunogen

Amino acids PDPYRQLREK were used as the immunogen for this TMEM205 antibody. **Storage** Aliquot and store the TMEM205 antibody at -20oC.