


TNK1 Antibody [clone 170CT7.5.4] (F40236)

Catalog No.	Formulation	Size
F40236-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F40236-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Mouse
Format	Purified
Host	Mouse
Clonality	Monoclonal (mouse origin)
Isotype	Mouse IgG1, k
Clone Name	170CT7.5.4
Purity	Purified
UniProt	Q13470
Applications	Western Blot : 1:100-1:1000
Limitations	This TNK1 antibody is available for research use only.



Description

TNK1 is involved in negative regulation of cell growth. Has tumor suppressor properties. Plays a negative regulatory role in the Ras-MAPK pathway. May function in signaling pathways utilized broadly during fetal development and more selectively in adult tissues and in cells of the lymphohematopoietic system. Could specifically be involved in phospholipid signal transduction.

Application Notes

Titration of the TNK1 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

This TNK1 antibody was produced from a mouse immunized with TNK1 recombinant protein.

Storage

Aliquot the TNK1 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.