


TRPV5 Antibody (R35921)

| Catalog No. | Formulation | Size |
|--------------|--|--------|
| R35921-100UG | 0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide $$ | 100 ug |

Bulk quote request

| Availability | 1-3 business days |
|----------------------|---|
| Species Reactivity | Human |
| Predicted Reactivity | Cow, Dog, Mouse, Pig, Rat |
| Format | Antigen affinity purified |
| Clonality | Polyclonal (goat origin) |
| Isotype | Goat Ig |
| Purity | Antigen affinity |
| Gene ID | 56302 |
| Applications | Western Blot : 1-3ug/ml ICC/IF : suitable ELISA (peptide) LOD : 1:16000 |
| Limitations | This TRPV5 antibody is available for research use only. |



Description

Additional name(s) for this target protein: Transient receptor potential cation channel, subfamily V, member 5

Application Notes

Optimal dilution of the TRPV5 antibody should be determined by the researcher.

Immunogen

Amino acids SHRGWEILRQNT were used as the immunogen for this TRPV5 antibody.

Storage

Aliquot and store the TRPV5 antibody at -20oC.