


TSPAN5 Antibody (R35904)

Catalog No.	Formulation	Size
R35904-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide $$	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Predicted Reactivity	Cow, Dog, Pig
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	10098
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This TSPAN5 antibody is available for research use only.



Description

Additional name(s) for this target protein: Tetraspanin 5, NET-4

Application Notes

Optimal dilution of the TSPAN5 antibody should be determined by the researcher.

Immunogen

 $\label{thm:local_signal_signal} Amino\ acids\ SGKHYKGPEVSC\ were\ used\ as\ the\ immunogen\ for\ this\ TSPAN5\ antibody.$ **Storage** Aliquot and store the TSPAN5 antibody at -20oC.