


TSPY3 Antibody (F40043)

Catalog No.	Formulation	Size
F40043-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F40043-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	P0CV98
Applications	Western Blot : 1:1000
Limitations	This TSPY3 antibody is available for research use only.



Description

The function of TSPY3 is unknown. TSPY3 may be involved in sperm differentiation and proliferation

Application Notes

Titration of the TSPY3 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 277-305 from the human protein was used as the immunogen for this TSPY3 antibody.

Storage

Aliquot the TSPY3 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.