


TUBD1 Antibody (Tubulin delta 1) (F52159)

Catalog No.	Formulation	Size
F52159-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F52159-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Primate
Format	Antigen affinity purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	Q9UJT1
Applications	Western Blot : 1:1000
Limitations	This TUBD1 antibody is available for research use only.



Description

In the elongating spermatid it is associated with the manchette, a specialized microtubule system present during reshaping of the sperm head (By similarity).

Application Notes

Titration of the TUBD1 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 366-392 from the human protein was used as the immunogen for this TUBD1 antibody.

Storage

Aliquot the TUBD1 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.