


## U2AF1L4 Antibody (R33731)

Catalog No.	Formulation	Size
R33731-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	199746
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This U2AF1L4 antibody is available for research use only.



Western blot testing of human cerebellum lysate with U2AF1L4 antibody at 0.1ug/ml.  
Predicted molecular weight: ~21kDa.

## Description

Additional name(s) for this target protein: U2 small nuclear RNA auxiliary factor 1-like 4

## Application Notes

Optimal dilution of the U2AF1L4 antibody should be determined by the researcher.

## Immunogen

Amino acids EDGERAVAELSNR were used as the immunogen for this U2AF1L4 antibody.

## **Storage**

Aliquot and store the U2AF1L4 antibody at -20oC.