


UCHL5 Antibody (F47733)

Catalog No.	Formulation	Size
F47733-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F47733-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	Q9Y5K5
Applications	Western Blot : 1:1000
Limitations	This UCHL5 antibody is available for research use only.



Description

Protease that specifically cleaves 'Lys-48'-linked polyubiquitin chains. Deubiquitinating enzyme associated with the 19S regulatory subunit of the 26S proteasome.

Application Notes

Titration of the UCHL5 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 136-164 from the human protein was used as the immunogen for this UCHL5 antibody.

Storage

Aliquot the UCHL5 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.