


# UNC5B Antibody (R35840)

Catalog No.	Formulation	Size
R35840-100UG	$0.5\ \text{mg/ml}$ in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

#### **Bulk quote request**

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Pig, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	219699
Applications	Western Blot : 2-4ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This UNC5B antibody is available for research use only.



Western blot testing of human bone marrow lysate with UNC5B antibody at 2ug/ml with [B] and without [A] blocking/immunizing peptide. Predicted molecular weight: 104kDa.

### **Description**

Additional name(s) for this target protein: Unc-5 homolog B

## **Application Notes**

Optimal dilution of the UNC5B antibody should be determined by the researcher.

#### **Immunogen**

Amino acids SARDWIFQLKTQAHQ were used as the immunogen for this UNC5B antibody. **Storage** Aliquot and store the UNC5B antibody at -20oC.