

Urokinase Antibody (F52084)

Catalog No.	Formulation	Size
F52084-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F52084-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Format	Antigen affinity purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	P06869
Applications	Western Blot : 1:1000
Limitations	This Urokinase antibody is available for research use only.

95
72
55
36
28

Urokinase antibody western blot analysis in mouse lung tissue lysate.

95
72
55
36
28

Urokinase antibody western blot analysis in human placenta tissue lysate

Description

Specifically cleave the zymogen plasminogen to form the active enzyme plasmin.

Application Notes

Titration of the Urokinase antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 304-331 from the mouse protein was used as the immunogen for this Urokinase antibody.

Storage

Aliquot the Urokinase antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.