


USP32 Antibody (F53864)

Catalog No.	Formulation	Size
F53864-0.2ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.2 ml
F53864-0.05ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.05 ml

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	Q8NFA0
Applications	Western Blot : 1:500-1000
Limitations	This USP32 antibody is available for research use only.



Description

USP32 (ubiquitin specific peptidase 32) catalyzes the conversion of a ubiquitin C-terminal thioester to a free ubiquitin and a thiol. It contains one DUSP domain and three EF-hand calcium binding domains and is localized to membranes in a lipid-anchored fashion.

Application Notes

Titration of the USP32 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 438-470 from the human protein was used as the immunogen for the USP32 antibody. **Storage** Aliquot the USP32 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles. Ordering: Phone: 858.663.9055 | Fax: 1.267.821.0800 | Email: info@nsjbio.com Copyright $\ensuremath{\texttt{@}}$ NSJ Bioreagents. All rights reserved.