


Zebrafish Thrb Antibody (R33481)

Catalog No.	Formulation	Size
R33481-100UG	0.5~mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Zebrafish
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Applications	Western Blot : 2-5ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This zebrafish Thrb antibody is available for research use only.



Western blot testing of three different zebrafish adult heart lysate samples with Thrb antibody at 2.5ug/ml. Predicted molecular weight: ~44kDa.

Description

Additional name(s) for this target protein: Thyroid hormone receptor beta

Application Notes

Optimal dilution of the Thrb antibody should be determined by the researcher.

Immunogen

Amino acids DIGSAPIVNAPEGNKVD were used as the immunogen for this zebrafish Thrb antibody.

Storage

