


# ZIC2 Antibody (R34286)

Catalog No.	Formulation	Size
R34286-100UG	0.5~mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

#### **Bulk quote request**

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	7546
Applications	Western Blot: 0.5-1.5ug/ml ELISA (peptide) LOD: 1:1000
Limitations	This ZIC2 antibody is available for research use only.



RWPE1 cell lysate overexpressing human ZIC2-MYC tag tested with ZIC2 antibody (0.5ug/ml) in the right panel and probed with anti-MYC tag (1/1000) and anti-beta-Actin in the left panel. Mock-transfected RWPE1 lysate in lane 1 and expressing GFP in lanes

### **Description**

Additional name(s) for this target protein: Zic family member 2

## **Application Notes**

Optimal dilution of the ZIC2 antibody should be determined by the researcher.

#### **Immunogen**

Amino acids HSGLSSNFNEWY were used as the immunogen for this ZIC2 antibody. **Storage** Aliquot and store the ZIC2 antibody at -20oC.