


## ABHD8 Antibody (F41117)

Catalog No.	Formulation	Size
F41117-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F41117-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

<b>Availability</b>	1-3 business days
<b>Species Reactivity</b>	Human
<b>Predicted Reactivity</b>	Primate
<b>Format</b>	Antigen affinity purified
<b>Clonality</b>	Polyclonal (rabbit origin)
<b>Isotype</b>	Rabbit Ig
<b>Purity</b>	Antigen affinity
<b>UniProt</b>	Q96113
<b>Applications</b>	Western blot : 1:1000
<b>Limitations</b>	This ABHD8 antibody is available for research use only.



## Description

This gene is upstream of, and in a head-to-head orientation with the gene for the mitochondrial ribosomal protein L34. Abhydrolase domain-containing protein 8 contains alpha/beta hydrolase fold and secretory lipase domains.

## Application Notes

Titration of the ABHD8 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

## Immunogen

A portion of amino acids 410-439 from the human protein was used as the immunogen for this ABHD8 antibody.

## Storage

Aliquot the ABHD8 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.