


TLE2 Antibody (F42039)

Catalog No.	Formulation	Size
F42039-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F42039-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Format	Antigen affinity purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	Q04725
Applications	Western Blot : 1:1000
Limitations	This TLE2 antibody is available for research use only.



Description

TLE2 is a transcriptional corepressor that binds to a number of transcription factors. Inhibits the transcriptional activation mediated by CTNNB1 and TCF family members in Wnt signaling. The effects of full-length TLE family members may be modulated by association with dominant-negative AES (By similarity).

Application Notes

Titration of the TLE2 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 637-665 from the human protein was used as the immunogen for this TLE2 antibody.

Storage

Aliquot the TLE2 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.