


NKX3.2 Antibody (F45860)

Catalog No.	Formulation	Size
F45860-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F45860-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Xenopus
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	P78367
Applications	Western blot : 1:1000
Limitations	This NKX3.2 antibody is available for research use only.



Description

This gene encodes a member of the NK family of homeobox-containing proteins. The encoded protein may play a role in skeletal development.

Application Notes

Titration of the NKX3.2 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 3-2 from the human protein was used as the immunogen for this NKX3.2 antibody.

Storage

Aliquot the NKX3.2 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.