


## AATK Antibody (F53414)

Catalog No.	Formulation	Size
F53414-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F53414-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

<b>Availability</b>	1-3 business days
<b>Species Reactivity</b>	Human, Mouse
<b>Format</b>	Antigen affinity purified
<b>Clonality</b>	Polyclonal (rabbit origin)
<b>Isotype</b>	Rabbit Ig
<b>Purity</b>	Antigen affinity
<b>UniProt</b>	Q6ZMQ8
<b>Applications</b>	Western blot : 1:2000
<b>Limitations</b>	This AATK antibody is available for research use only.



## Description

Apoptosis-associated tyrosine kinase may be involved in neuronal differentiation.

## Application Notes

Titration of the AATK antibody may be required due to differences in protocols and secondary/substrate sensitivity.

## Immunogen

This AATK antibody was produced from a rabbit immunized with a KLH conjugated synthetic peptide between 55-90 amino acids from the human AATK.

## Storage

Aliquot the AATK antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.