

NDUFS8 Antibody / NADH-ubiquinone oxidoreductase 23 kDa subunit (R33125)

Catalog No.	Formulation	Size
R33125-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Rat, Dog, Cow
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	4728
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:2000
Limitations	This NDUFS8 antibody is available for research use only.



Western blot testing of mouse heart lysate with NDUFS8 antibody at 0.1ug/ml. Predicted molecular weight: ~24kDa.

Description

Additional name(s) for this target protein: NADH dehydrogenase (ubiquinone) Fe-S protein 8

Application Notes

Optimal dilution of the NDUFS8 antibody should be determined by the researcher.

Immunogen

Amino acids AEPRADGSRR were used as the immunogen for this NDUFS8 antibody.

Storage

Aliquot and store the NDUFS8 antibody at -20oC.