


Acylglycerol kinase Antibody (R33214)

Catalog No.	Formulation	Size
R33214-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	55750
Applications	Western blot : 1-3ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This Acylglycerol kinase antibody is available for research use only.



Description

Additional name(s) for this target protein: AGK; MULK; multi-substrate lipid kinase; AGK

Application Notes

Optimal dilution of the Acylglycerol kinase antibody should be determined by the researcher.

Immunogen

Amino acids CDPRKREQMLTSP were used as the immunogen for this Acylglycerol kinase antibody.

Storage

Aliquot and store the Acylglycerol kinase antibody at -20oC.