


## Agbl5 Antibody [isoform 3] (R33232)

Catalog No.	Formulation	Size
R33232-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

<b>Availability</b>	1-3 business days
<b>Species Reactivity</b>	Mouse
<b>Format</b>	Antigen affinity purified
<b>Clonality</b>	Polyclonal (goat origin)
<b>Isotype</b>	Goat Ig
<b>Purity</b>	Antigen affinity
<b>Gene ID</b>	231093
<b>Applications</b>	Western blot : 1-3ug/ml ELISA (peptide) LOD : 1:16000
<b>Limitations</b>	This Agbl5 antibody is available for research use only.



## Description

Additional name(s) for this target protein: ATP/GTP binding protein-like 5

## Application Notes

Optimal dilution of the Agbl5 antibody should be determined by the researcher.

1. This Agbl5 antibody is specific to isoform 3 only.

## Immunogen

Amino acids CFSKPEEAGSHVE were used as the immunogen for this Agbl5 antibody. (1)

## Storage

Aliquot and store the Agbl5 antibody at -20oC.