


AIM2 Antibody (R33510)

Catalog No.	Formulation	Size
R33510-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	9447
Localization	Cytoplasmic, membrane, nuclear
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This AIM2 antibody is available for research use only.



Western blot testing of Daudi lysate with AIM2 antibody at 0.3ug/ml. Predicted molecular weight ~39kDa.

Description

Additional name(s) for this target protein: Absent in melanoma 2

Application Notes

Optimal dilution of the AIM2 antibody should be determined by the researcher.

Immunogen

Amino acids DKQYKSVTKPKPLSQ were used as the immunogen for this AIM2 antibody.

Storage

Aliquot and store the AIM2 antibody at -20oC.