


## ABCB9 Antibody (R34095)

Catalog No.	Formulation	Size
R34095-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

<b>Availability</b>	1-3 business days
<b>Species Reactivity</b>	Human
<b>Format</b>	Antigen affinity purified
<b>Clonality</b>	Polyclonal (goat origin)
<b>Isotype</b>	Goat Ig
<b>Purity</b>	Antigen affinity
<b>Gene ID</b>	23457
<b>Applications</b>	Immunohistochemistry (FFPE) : 3-5ug/ml ELISA (peptide) LOD : 1:64000
<b>Limitations</b>	This ABCB9 antibody is available for research use only.



IHC testing of FFPE human testis with ABCB9 antibody at 4ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

## Description

Additional name(s) for this target protein: ATP-binding cassette, sub-family B (MDR/TAP), member 9

## Application Notes

Optimal dilution of the ABCB9 antibody should be determined by the researcher.

## Immunogen

Amino acids GHNEPVANGSHKA were used as the immunogen for this ABCB9 antibody.

## Storage

Aliquot and store the ABCB9 antibody at -20oC.