


## AMHR2 Antibody (R34253)

| Catalog No.  | Formulation   | Size   |
|--------------|---|--------|
| R34253-100UG | 0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide | 100 ug |

[Bulk quote request](#)

|                           |   |
|---------------------------|---|
| <b>Availability</b>       | 1-3 business days   |
| <b>Species Reactivity</b> | Human   |
| <b>Format</b>             | Antigen affinity purified                                 |
| <b>Clonality</b>          | Polyclonal (goat origin)                                  |
| <b>Isotype</b>            | Goat Ig   |
| <b>Purity</b>             | Antigen affinity  |
| <b>Gene ID</b>            | 269   |
| <b>Applications</b>       | Western blot : 1-3ug/ml<br>ELISA (peptide) LOD : 1:128000 |
| <b>Limitations</b>        | This AMHR2 antibody is available for research use only.   |



## Description

Additional name(s) for this target protein: AMH type II receptor, Mullerian hormone receptor, type II

## Application Notes

Optimal dilution of the AMHR2 antibody should be determined by the researcher.

1. This AMHR2 antibody is specific for isoforms 1 and 3.

## Immunogen

Amino acids HPQESHFPESCPR were used as the immunogen for this AMHR2 antibody. (1)

## Storage

Aliquot and store the AMHR2 antibody at -20oC.