


MIA Antibody (R35042)

Catalog No.	Formulation	Size
R35042-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	8190
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This MIA antibody is available for research use only.



Western blot testing of human glioma lysate with MIA antibody at 1ug/ml. Predicted molecular weight: ~14kDa.

Description

Additional name(s) for this target protein: Melanoma inhibitory activity, CD-RAP

Application Notes

Optimal dilution of the MIA antibody should be determined by the researcher.

Immunogen

Amino acids PKLADRKLCADQE were used as the immunogen for this MIA antibody.

Storage

Aliquot and store the MIA antibody at -20oC.