



## ARHGEF4 Antibody (R35175)

Catalog No.	Formulation	Size
R35175-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

<b>Availability</b>	1-3 business days
<b>Species Reactivity</b>	Human
<b>Format</b>	Antigen affinity purified
<b>Clonality</b>	Polyclonal (goat origin)
<b>Isotype</b>	Goat Ig
<b>Purity</b>	Antigen affinity
<b>Gene ID</b>	50649
<b>Applications</b>	Western blot : 2-4ug/ml ELISA (peptide) LOD : 1:32000
<b>Limitations</b>	This ARHGEF4 antibody is available for research use only.





## Description

Additional name(s) for this target protein: APC-stimulated guanine nucleotide exchange factor 4

## Application Notes

Optimal dilution of the ARHGEF4 antibody should be determined by the researcher.

## Immunogen

Amino acids PWEEPAGEKPCS were used as the immunogen for this ARHGEF4 antibody.

## Storage

Aliquot and store the ARHGEF4 antibody at -20°C.