



ASXL1 Antibody (RQ4183)

Catalog No.	Formulation	Size
RQ4183	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug


Bulk quote request

Availability	1-3 business days
Species Reactivity	Human, Mouse
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit IgG
Purity	Antigen affinity purified
Buffer	Lyophilized from 1X PBS with 2% Trehalose and 0.025% sodium azide
UniProt	Q8IXJ9
Applications	Western Blot : 0.5-1ug/ml Flow cytometry : 1-3ug/10 ⁶ cells
Limitations	This ASXL1 antibody is available for research use only.





Flow cytometry testing of human U-2 OS cells with ASXL1 antibody at $1\mu\text{g}/10^6$ cells (blocked with goat sera); Red=cells alone, Green=isotype control, Blue=ASXL1 antibody.



Flow cytometry testing of human HepG2 cells with ASXL1 antibody at $1\mu\text{g}/10^6$ cells (blocked with goat sera); Red=cells alone, Green=isotype control, Blue=ASXL1 antibody.

Description

Putative Polycomb group protein ASXL1 is a protein that in humans is encoded by the ASXL1 gene. This gene is similar to the *Drosophila* additional sex combs gene, which encodes a chromatin-binding protein required for normal determination of segment identity in the developing embryo. The protein is a member of the Polycomb group of proteins, which are necessary for the maintenance of stable repression of homeotic and other loci. The protein is thought to disrupt chromatin in localized areas, enhancing transcription of certain genes while repressing the transcription of other genes. The protein encoded by this gene functions as a ligand-dependent co-activator for retinoic acid receptor in cooperation with nuclear receptor coactivator 1. Mutations in this gene are associated with myelodysplastic syndromes and chronic myelomonocytic leukemia. Alternative splicing results in multiple transcript variants.

Application Notes

Optimal dilution of the ASXL1 antibody should be determined by the researcher.

Immunogen

Amino acids KKERTWAEARLVLENYSAPMTPKQILQVIEAE were used as the immunogen for the ASXL1 antibody.

Storage

After reconstitution, the ASXL1 antibody can be stored for up to one month at 4°C. For long-term, aliquot and store at -20°C. Avoid repeated freezing and thawing.