



APAF1 Antibody [clone BHF-1] (RQ4704)

Catalog No.	Formulation	Size
RQ4704	Antibody in PBS with 0.02% sodium azide and 50% glycerol	100 ug

[Bulk quote request](#)

Availability	1-2 weeks
Species Reactivity	Human, Mouse
Format	Purified
Clonality	Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	BHF-1
Purity	Affinity purified
UniProt	O14727
Applications	Western blot : 1:500-1:2000 Immunohistochemistry (FFPE) : 1:50-1:200
Limitations	This APAF1 antibody is available for research use only.





IHC staining of FFPE mouse kidney tissue with APAF1 antibody. HIER: boil tissue sections in pH6, 10mM citrate buffer, for 10-20 min and allow to cool before testing.

Description

Apoptotic Peptidase Activating Factor 1 is a cytoplasmic protein that initiates apoptosis. And it is an essential downstream effector of p53-mediated apoptosis.

Application Notes

Optimal dilution of the APAF1 antibody should be determined by the researcher.

Immunogen

A synthetic peptide specific to human Apoptotic Peptidase Activating Factor 1 was used as the immunogen for the APAF1 antibody.

Storage

Store the APAF1 antibody at -20oC.