


58K Golgi protein Antibody [clone AODH-6] (RQ4820)

Catalog No.	Formulation	Size
RQ4820	Antibody in PBS with 0.02% sodium azide, 50% glycerol and 0.4-0.5mg/ml BSA	100 ul

[Bulk quote request](#)

Availability	1-2 weeks
Species Reactivity	Human
Format	Purified
Host	Rabbit
Clonality	Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	AODH-6
Purity	Affinity purified
UniProt	O95954
Applications	Western Blot : 1:500-1:2000
Limitations	This 58K Golgi protein antibody is available for research use only.



Western blot testing of human fetal liver with 58K Golgi protein antibody. Predicted molecular weight ~58 kDa.

Description

Formimidoyltransferase-cyclodeaminase (FTCD), also called '58K Golgi protein' binds and promotes bundling of vimentin filaments originating from the Golgi.

Application Notes

Optimal dilution of the 58K Golgi protein antibody should be determined by the researcher.

Immunogen

A synthetic peptide specific to human 58K Golgi protein / FTCD was used as the immunogen for the 58K Golgi protein antibody.

Storage

Store the 58K Golgi protein antibody at -20oC.