


## 58K Golgi protein Antibody / Golgi Marker (R36023)


Catalog No.	Formulation	Size
R36023-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Pig
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	10841
Applications	Western Blot : 0.03-0.5ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This 58K Golgi protein antibody is available for research use only.



Western blot testing of human liver lysate with 58K Golgi protein antibody at 0.03ug/ml.  
Predicted molecular weight: ~58 kDa.



Western blot testing of human liver lysate with 58K Golgi protein antibody at 0.05ug/ml.  
Predicted molecular weight: ~58 kDa.

## Description

Additional name(s) for this target protein: FTCD, formiminotransferase cyclodeaminase

## Application Notes

Optimal dilution of the 58K Golgi protein antibody should be determined by the researcher.

## Immunogen

Amino acids SQLVECVPNFSEGKNQ were used as the immunogen for this 58K Golgi protein antibody.

## Storage

Aliquot and store the 58K Golgi protein antibody at -20oC.