


## 58K Golgi protein Antibody / Golgi Marker (R36023)


Catalog No.	Formulation	Size
R36023-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Pig
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	10841
Applications	Western Blot : 0.03-0.5ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This 58K Golgi protein antibody is available for research use only.



Western blot testing of human liver lysate with 58K Golgi protein antibody at 0.03ug/ml.  
Predicted molecular weight: ~58 kDa.



Western blot testing of human liver lysate with 58K Golgi protein antibody at 0.05ug/ml.  
Predicted molecular weight: ~58 kDa.

## Description

Additional name(s) for this target protein: FTCD, formiminotransferase cyclodeaminase

## Application Notes

Optimal dilution of the 58K Golgi protein antibody should be determined by the researcher.

## Immunogen

Amino acids SQLVECVPNFSEGKNQ were used as the immunogen for this 58K Golgi protein antibody.

## Storage

Aliquot and store the 58K Golgi protein antibody at -20oC.