


## 5HT3B receptor Antibody (R35996)

Catalog No.	Formulation	Size
R35996-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	9177
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This 5HT3B receptor antibody is available for research use only.



Western blot testing of human cerebral cortex lysate with 5HT3B receptor antibody at 1ug/ml. Predicted molecular weight (human): ~50kDa.

### Description

Additional name(s) for this target protein: 5-hydroxytryptamine (serotonin) receptor 3B; HTR3B

### Application Notes

Optimal dilution of the Serotonin Receptor 3B / HTR3B antibody should be determined by the researcher.

## Immunogen

Amino acids SNYLQTQDQTDQQE were used as the immunogen for this 5HT3B receptor antibody.

## Storage

Aliquot and store the 5HT3B receptor antibody at -20oC.